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SOFTWARE COMPATIBILITY

Mightex	PolyScan2 software included in every package free of charge
3rd Party Support	Nikon's NIS Elements (Polygon400-G model only)
	NeuroPG: Open Source Software for Optical Pattern Generation
	Micro-Manager Open Source Microscopy Software

SYSTEM AND COMMUNICATION


Operating System	Windows XP, Vista, 7, 8, and 10
Interface	USB2.0
Power Supply	5Vdc 3A input power
Screen Resolution	1,366x768 or higher

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
<https://www.mightexbio.com/products/polygon400/polygon400-2/>


Photostimulation: Neuroscience Optogenetics

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
Photostimulation: Cell Biology Optogenetics

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



Photostimulation: Photo-**activation**, conversion, switching

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Photostimulation: Photo-activation, **conversion**, switching

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Photostimulation: Photo-activation, conversion, **switching**

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





FRAP
(Fluorescence Recovery After Photobleaching)

iFRAP
(Inverse FRAP)

FLIP
(Fluorescence Lost in Photobleaching)

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



自由行動下の Optogenetics Deep-Brain, Cortex

Calcium Imaging + Single-Cell
Resolution Optogenetics

Calcium Imaging +
Widefield Optogenetics

Calcium Imaging

Calcium Imaging +
Targeted Optogenetics

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




Single-Cell Resolution Optogenetic Spiking

- (a) Illustration of highly targeted optical stimulation of a single ChR2-mCherry expressing neuron in the mouse somatosensory cortex. Using a small ($25\text{ }\mu\text{m}$ diameter), low-powered (3 mW) spot of illumination centred on the target cell, action potentials could be induced in current clamp (top trace), without any indication in voltage clamp of post-synaptic currents caused by spiking in other neurons (bottom trace).
- (b) Illustration of illumination with a large ($125\text{ }\mu\text{m}$ diameter), high-powered (15 mW) spot. Multiple spikes were induced in current clamp, and voltage clamp traces showed evidence of post-synaptic currents caused by spiking in other neurons.
- Courtesy **Matthew Tran & Dr Blake Richards**,
University of Toronto, Canada.


ChR2 Assisted Circuit Mapping (CRACM)

- (A) Image of an acute brain slice prepared from a Thy1-ChR2-EYFP mouse with ChR2 expression in L5 pyramidal neurons. Whole-cell patch clamp recording from a L2/3 cell.
(B) Light-activated excitatory postsynaptic potentials (EPSPs) triggered by patterned illumination of a 10×10 grid with 473 nm LED.
(C) Colormap of the activation pattern.
(D) Histogram of automatically measured responses from all cells in a grid. Objective lens is 10X.
- Courtesy **Qiuyu Wu & Dr Alexander Chubykin**, *Purdue University, USA*.


Optogenetic Stimulation of Synaptic Inputs

- Light-evoked inhibitory postsynaptic current (IPSC) recorded in a transgenic mice expressing ChR2 in GABAergic neurons. The postsynaptic cell is non-GABAergic (ChR2 negative) and blue light stimulates GABAergic afferents expressing ChR2. Blue bars indicate the time of light illumination. Spot 1 illuminated by Polygon400 evoked reliable IPSCs whereas Spot 2 caused no response.
- Courtesy **Dr Wataru Inoue**,
University of Western Ontario,
Canada.


Mapping Optically Induced Depolarization in ChR2-expressing Hippocampal Neurons using Polygon400

- E18 Sprague-Dawley rat neurons were transduced with CamKII-ChR2-GFP lentivirus. Somatic activity was recorded via whole-cell patch-clamp electrophysiology. Each field was illuminated by the Polygon400 at 100% power and 20ms exposure time. Intensity of magenta pattern represents depolarization with respect to the instantaneous resting potential prior to stimulation with the Polygon's 470nm LED. In the image, green represents the magnitude of GFP signal and black represents the fluorescence intensity of AlexaFluor 594 backfilled by the patch pipette.
- Courtesy **Dr Jacob Robinson**, *Rice University, USA*.


Photopatterning of an Acrylate Film

- A 50 μm film of liquid tetra (ethylen glycol) diacrylate, containing 1 wt% Irgacure 819 photoinitiator was irradiated through a glass coverslip, using the Polygon400 and a 400nm LED light source to project a pattern onto the film surface. (A) Projection of CU logo image using the 400nm LED and a 4X objective. (B) Brightfield image of the resulting pattern in the film. Photopolymerization causes a large change in refractive index in the resin, allowing immediate visualization of the pattern. Standard development techniques could subsequently be performed by washing the film in solvent to remove the unexposed areas of liquid resin.
- Courtesy **Gayla Berg**, *University of Colorado, USA*.


Achieving Subcellular Resolution of Input Activation in a 250um-thick Acute Cortical Section Using the Polygon400

- NMDA-receptor mediated synaptic currents recorded at a +40mV membrane potential were elicited by light activation of channelrhodopsin-expressing terminals of thalamocortical afferents onto an upper layer cortical GABAergic interneuron. moving the 470nm rectangular illumination to the right or the left of the cell activates inputs of different strength innervating different somatodendritic domains of the cell. In black and red are the individual and the average traces respectively.
- Courtesy **Dr. Theofanis Karayannis**, *University of Zurich, Switzerland*.


In Vivo Optogenetic Control of Zebrafish Larva Using Polygon400 Illuminator

- A) Light-evoked response of a head-fixed larva expressing channelrhodopsin (right). Photostimulation site was indicated by a blue circle (470nm).
- B) A schematic diagram of the Polygon400 patterned illumination in in vivo optogenetic mapping system.
- Courtesy **Dr. Sachiko Tsuda**, *Saitama University*, Japan.


Uncaging RuBi Glutamate with the Polygon400 to Study Neuronal Networks in the Cerebellar Cortex

- An example of connectivity mapping that allow to reproduce some results from in Valera et al. (eLife, 2016). 100µm RuBi glutamate was uncaged at various locations in the granular layer with 20ms pulses (blue bar), exciting notably granule cells. We can then measure the spatial organization of the granule cells to Purkinje cell (PC) connections by recording PCs in whole cell patch clamp. Granule cells triggers both monosynaptic excitatory current onto PCs (left map, measured at -60mV), but also disynaptic inhibitory currents via molecular layer interneurons (right map, measured at 0 mV). Average evoked responses at one location (dotted blue square) are shown at the bottom of the figure. Data storage, measurements, and map representations can be made online, using a homemade software in python.
- Courtesy, **Dr. Antoine Valera** & **Dr. Angus Silver**, University College London, UK.


Optogenetic Activation of Channelrhodopsin in Transfected Hippocampal Neurons Using Mightex's Polygon400

- A) Patch-clamp recording of the current generated by the channelrhodopsin when activated with the blue LED light from the Polygon400 illuminator as indicated by the blue bar.
- B) Activation of an action potential (upper trace) and the corresponding current under voltage-clamp conditions (lower trace). The action potential could be evoked by a 0.5ms stimulation of the blue LED light from the Polygon400 illuminator at 60% intensity.
- Courtesy, **Dr. Hans van Hooft**, *University of Amsterdam*, Netherlands.


ing Polygon400 to Stimulate Hippocampal Slices Acutely Prepared from Mice Transgenically Expressing ChR2 in the Dentate Gyrus

- The first figure a repetitive stimulation with a larger circular pattern which elicited increasing responses until an action potential is generated. The second figure is increasing the intensity of stimulation from 30% to 100% and keeping the size of the pattern the same. Here we're driving the cell at 10Hz at 100% and one can see that the stimulation is sufficient to drive doublets of action potentials during each bout of depolarization.
- Courtesy, **Dr. Geoffrey G. Murphy**,
University of Michigan, USA


Local Photostimulation of Channelrhodopsin-2 Using Polygon400 Illuminator

- A) Acute brain slice from a YFP-channelrhodopsin-2 (ChR-2) mouse depicting its expression in cortical L5 pyramidal neurons. B) L5 pyramidal neuron from the somatosensory cortex filled with Alexa 594 to allow visualization of neuronal compartments without stimulating ChR-2. Blue dots indicate on the illuminated areas (Blue LED – 470 nm) along the apical dendrite. Expansion of the marked areas depicting the delicate dendrites that were stimulated. C) Electrophysiological (patch clamp) current recordings from the soma, corresponding to local photostimulation of ChR-2 by blue light. The numbers in the bottom of the trace corresponds to the stimulated dot numbers as indicated in B.
- Courtesy, **Dr. Yossi Buskila**, *University of Western Sydney, Australia.*

Purkinje Cell Firing Properties Unveiled by Optogenetic Activation of the Cerebellar Granular Layer Using Variable Light Patterns with Mightex Polygon400

- Purkinje cell (PC) firing is monitored, while various light stimulation patterns are delivered to the granular layer (gl) and/or molecular layer (ml) of acute mouse cerebellar slices. The mice cerebellum was injected with AAVs carrying the ChR2 construct. The figure shows the firing frequency increase caused in a PC unit (scale bar 500ms) by optogenetic activation of a granular layer ROI (blue rectangle).
- Courtesy, **Dr. Lisa Mapelli & Dr. Simona Tritto**, *University of Pavia, Italy*.

Polygon400 - 文献 2018

- [Che, A., Babij, R., Iannone, A.F., Fetcho, R.N., Ferrer, M., Liston, C., Fishell, G., & De Marco Garcia, N.V. "Layer I Interneurons Sharpen Sensory Maps during Neonatal Development". **Neuron**, \(2018\).](#)
- [Gomati, S.V., Schafer, C.B., Eelkman Rooda, O.H.J., Nigg, A.L., De Zeeuw, C.I., & Hoebeek, F.E. \(2018\). "Differentiating Cerebellar Impact on Thalamic Nuclei". **Cell Reports**, 23\(9\), 2690-2704, \(2018\)](#)
- [Dine, E., Gil, A.A., Uribe, G., Brangwyne C.P., & Toettcher J.E. "Protein Phase Separation Provides Long-Term Memory Transient Spatial Stimuli". **Cell Systems**, \(2018\).](#)
- [Oboti, L., Russo, E., Tran, T., Durstewitz, D., & Corbin, J.G. "Amygdala Corticofugal Input Shapes Mitral Cell Responses in Accessory Olfactory Bulb". **eNeuro**, \(2018\).](#)
- [van Haren, J., Charafeddine, R.A., Ettinger, A., Wang, H., & Wittmann, T. "Local Control of Intracellular Microtubule Dynamics by EB1 Photodissociation". **Nature Cell Biology**, 20, 252-261 \(2018\).](#)
- [Butler, J.L, Hay, Y.A., & Paulsen O. "Comparison of Three Gamma Oscillations in the Mouse Entorhinal-Hippocampal System". **European Journal of Neuroscience**, \(2018\).](#)

Polygon400 - 文献 2017

- [Andrasi, T., Verses, J.M., Rovira-Esteban, L., Kozma, R., Viktor, A., Gregori, E., & Hajos, N. "Differential excitatory control of 2 parallel basket cell networks in amygdala microcircuits". **PLoS Biology**, 15\(5\) \(2017\).](#)
- [Shimizu, K. & Stopfer, M. "A Population of Projection Neurons that Inhibits the Lateral Horn but Excites the Antennal Lobe through Chemical Synapses in Drosophila". **Frontiers in Neural Circuits**, 3\(11\) \(2017\).](#)
- [Adrian, M., Nijenhuis, W., Hoogstraaten, R.I., Willems, J., & Kapitein, L.C. "A Phytochrome-Derived Photoswitch for Intracellular Transport". **ACS Synthetic Biology**, 1248-1256 \(2017\).](#)
- [McBride, M.K., Hendriks, M., Liu, D.Q., Worrell, B.T., Broer, D.J., & Bowman, C.N. "Photoinduced Plasticity in Cross-Linked Liquid Crystalline Networks". **Advanced Materials**, 29\(17\) \(2017\).](#)
- [Malyshev, A.Y., Roshchin, M.V., Smirnova, G.R., Dolgikh, D.A., Balaban, P.M., & Ostrovsky, M.A. "Chloride conducting light activated channel GtACR2 can produce both cessation of firing and generation of action potentials in cortical neurons in response to light". **Neuroscience Letters**, 640, 76-80 \(2017\).](#)
- [Johnson, H.E., Goyal, Y., Pannucci, N.L., Schupbach, T., Shvarstman, S.Y., & Toettcher, J.E. "The spatiotemporal limits of Erk signalling". **Developmental Cell**, 40\(2\), 185-192 \(2017\).](#)
- [Duan, J.G., Fu, H., & Zhang, J.Y. "Activation of Parvalbumin-Positive Neurons in Both Retina and Primary Visual Cortex Improves the Feature-Selectivity of Primary Visual Cortex Neurons". **Neuroscience Bulletin**, 33\(3\), 255-263 \(2017\).](#)
- [Watanabe, M., Feola, I., Majumder, R., Jangsangthong, W., Teplenin, A.S., Ypey, D.L., Schali, M.J., Zeppenfeld, K., Vries, A.A.F., & Pijnappels, D.A. "Optogenetic manipulation of anatomical re-entry by light-guided generation of a reversible local conduction block". **Cardiovascular Research**, 113\(3\), 354-366 \(2017\).](#)

Polygon400 - 文献 2014-2016

- [Fiedler, C.I., Aisenbrey, E.A., Heveran, C.M., Ferguson, V.L., Bryant, S.J., & McLeod, R.R. "Enhanced mechanical properties of photo-clickable thiol-ene PEG hydrogels through repeated photopolymerization of in-swollen macromer". **Soft Matter**, 44 \(2016\).](#)
- [Butler, J.L., Mendonca, P.R. F., Robinson, H.P.C., & Paulsen, O. "Intrinsic Corneal Ammonia Area 1 Theta-Nested Gamma Oscillations Induced by Optogenetic Theta Frequency Stimulation". **Journal of Neuroscience**, 36\(15\), 4155-4169 \(2016\).](#)
- [Konetski, D., Gong, T., & Bowman, C.N. "Photoinduced Vesicle Formation via the Copper-Catalyzed Azide-Alkyne Cycloaddition Reaction". **Langmuir**, 8195-8201 \(2016\).](#)
- [García, N.V.D.M., Priya, R., Tuncdemir, S.N., Fishell, G., & Karayannis, T. "Sensory inputs control the integration of neurogliaform interneurons into cortical circuits". **Nature Neuroscience**, 18, 393-401 \(2015\).](#)
- [Avants, B.W., Murphy, D.B., Dapello, J.A., & Robinson, J.T. "NeuroPG: open source software for optical pattern generation and data acquisition". **Frontiers in Neuroengineering** \(2015\).](#)
- [Peng, H.Y., Wang, C., Xi, W.X., Kowalski, B.A., Gong, T., Xie, X.L., Wang, W.T., Nair, D.P., McLeod, R.R., & Bowman, C.N. "Facile Image Patterning via Sequential Thiol-Michael/Thiol-Yne Click Reactions". **Chemistry of Materials**, 6819-6826 \(2014\).](#)